



Joson-Care 2024




Safety Quality Sustainability






Business Concept




History



Established in 2000, Joson-Care Co., Ltd. is Taiwan's premier manufacturer of hospital beds and ward equipment. We uphold rigorous quality control standards to deliver exceptional accuracy, durability, and reliability, while maximizing the use of renewable resources.





So far, we have sold our products to more than 40 countries. Such as U.S.A, Japan, France, South Korea, India, Iran, Panama, Thailand, Turkey and U.A.E.

Joson-Care's manufacturing system is meticulously designed to adhere to ISO 13485 and Quality Management System (QMS) standards. This rigorous approach ensures that our processes consistently deliver high-quality products, precisely meeting and often exceeding customer expectations for safety, performance, and reliability.



ICU Bed



19 series

- ES-19HD ICU Bed 7
- ES-19DF ICU Bed 8



12 series

- ES-12DW ICU Bed with Weight Scale 13
- ES-12HD Electric Hospital Bed 15
- ES-12DF Electric Hospital Bed 16
- ES-11DF Electric Hospital Bed 17



09.96 series

- ES-09DW Electric Hospital Bed with Weight Scale 21
- ES-96HD Electric Hospital Bed 21

Hospital Bed



- ES-05HD Electric Hospital Bed 25
- ES-05HDS Electric Hospital Bed 26
- MS-05HS Manual Hospital Bed 27
- ES-99HDS Electric Hospital Bed 27



02.16 series

- ES-02JH Electric Nursing Bed 31
- ES-16N series Electric Nursing Bed 32
- MS-16NS Manual Hospital Bed 33

Emergency Stretcher



Stretcher series

- JE-350 & JE-300 Hydraulic Stretcher 37
- JE-200 & JE-150 Hydraulic Stretcher 39
- JE-130 3 Manual Cranks Emergency Stretcher 42
- JE-120 2 Manual Cranks Emergency Stretcher 42
- JE-100 Simple Emergency Stretcher 42

43

Pediatric Bed



Pediatric Bed series

- JC-012 Electric Pediatric Bed 45
- JC-011 Manual Pediatric Bed 46
- JC-031 Infant Crib 47
- JC-030 Crib 47

Medical & Ward Equipment

Mattress	67
JT-100 Orthopedic Traction Frame	69
IV Stand	70
JA-107 Emergency Cart & UD Medical Cart	71
Medical Cart	73
Cabinet & Outpatient Service Equipment	74
Overbed Table	75
Caregiver Recliner & Doctor Chair	76

Home Care Bed

49



18 series

- ES-18UDS Ultra Low Bed 51
- ES-17UDS Ultra Low Bed 52



EN series

- EN-3M Foldable Electric Home Care Bed 57
- EN-2M Assembly Electric Home Care Bed 58




08 series

- ES-08FDS Electric Hospital Bed 63
- ES-08FSS Electric Hospital Bed(Metal frame) 64
- MS-08FS 3 Manual Cranks Hospital Bed 65
- MF-08FS 2 Manual Cranks Hospital Bed 65

19 series




ES-19HD / ES-19DF



Head Elevation
0 ~ 75°



Knee Elevation
0 ~ 30°



Bed Elevation



Auto-chair



Trend. $\pm 18^\circ$



Rev. Trend. $\pm 18^\circ$



Trend.
Angle Indicator




Head
Angle Indicator




Safety Standard



230KG
Work Load




4PCS
Actuator



Backup Battery



IP66
Actuator/Control
Box





Shear Stress
Bi-directional Release

ES-19HD

ICU Bed

Standards: IEC 60601-1, 60601-1-2, 60601-2-52



Shear Stress
Bi-directional Release



Bed Elevation 40cm



Backrest X-Ray
Module



Trend. / Rev. Trend. $\pm 18^\circ$



Head Angle Indicator /
Outer Control Panel





Night Light

ES-19DF

ICU Bed

Standards: IEC 60601-1, 60601-1-2, 60601-2-52



Chair Position



Bed Height Indicator



Brake Light Indicator



Storage Space



Bed-end Extension 20cm



Mattress Stopper


19 series

Shear Stress Bi-directional Release

The hip and back boards decline to minimize shear stress and ease pressure on the back and abdomen.



Bed Elevation 40cm



LO: 38 cm

● Green: In low pos.

HI: 78cm


● Orange: Above 45cm

Vertical Lifting

The bed's frame is fixed while elevation is being operate.



Non-moving Headboard



Headboard fixed during Trend, and elevation adjust.



Specification	ES-19HD	ES-19DF
Overall bed length	L : 2220mm	L : 2220mm
Overall bed width	W : 958mm	W : 958mm
(Siderails down/up)	990mm	990mm
Overall mattress base length	L : 1925mm	L : 1925mm
Overall mattress base width	W : 885mm	W : 885mm
Overall bed height	(5") H : 385~785mm (6") H : 410~810mm	(5") H : 385~785mm (6") H : 410~810mm
Head Elevation	0~75°	0~75°
Knee Elevation	0~30°	0~30°
Trend./Rev. Trend.	±18°	±18°
Work Load	230kg	230kg
Castor Diameter	5" or 6"	5" or 6"
Actuator	4	4
Volts & Frequency	AC 100-240V ; 50/60Hz	AC 100-240V ; 50/60Hz
Backup Battery	DC 24V 2.0Ah	DC 24V 2.0Ah
Actuator/Control Box	IP66	IP66
Safety Standard	IEC 60601-2-52	IEC 60601-2-52

Applicable mattress size (mm)	L1900xW890xH100	L1900xW890xH100
-------------------------------	-----------------	-----------------

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree














Features	ES-19HD	ES-19DF
Backrest X-Ray Module	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Caregiver Siderail Control	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Patient Siderail Control	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
LCD Screen Panel Remote Controller	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Backrest, Trend./Rev. Trend. Angle Indicator	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
CPR Quick Release	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Roller Bumper	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Lock-Out Function	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Auto Regression	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Examination Position	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Anti-Shock Position	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Chair Position	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Trend./Rev. Trend.	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Bed-end Extension	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Dust-Proof Central Braking/Steering Module	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
IV Pole Hole	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Drainage Bag Hanger	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Mattress Stopper	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Brake/Height Caution Indicator	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Night Light	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Third Siderail	<input type="checkbox"/>	<input type="checkbox"/>
Restraint Belt Hole	<input type="checkbox"/>	<input type="checkbox"/>

Standard Option














12 series



ES-12HD / ES-12DW

							
Head Elevation 0~75°	Knee Elevation 0~30°	Bed Elevation	Auto-chair	Trend. ±12°	Rev. Trend. ±12°	Trend. Angle Indicator	Head Angle Indicator
							
Safety Standard	Work Load 230KG	Actuator 4PCS	Backup Battery	Actuator/Control Box IP66			



ES-12DF / ES-11DF

							
Head Elevation 0~75°	Knee Elevation 0~30°	Bed Elevation	Auto-chair	Trend. ±12° (Option)	Rev. Trend. ±12° (Option)	Trend. Angle Indicator	Head Angle Indicator
							
Safety Standard	Work Load 230KG	Actuator 3PCS (Option)	Backup Battery	Actuator/Control Box IP66			

ES-12DW

ICU Bed
with Weight Scale

Standards: IEC 60601-1, IEC 60601-1-2, IEC 60601-2-52



Weight Monitoring Module



Bed Height Indicator



Trend. / Rev. Trend. $\pm 12^\circ$



Chair Position




Mattress Stopper

11 / 12 series


Safety Standard

11 & 12 series comply to latest safety standard IEC 60601-2-52, which ensures that all gaps between head/foot boards and siderails are smaller than 60mm



ES-12DW

Weight Scale & ACP Functions



Fall Prevention Alarm





Weight Sensors

Sensors placed around the bed to detect when a patient is attempting to get out of bed

ES-12HD

Electric Hospital Bed

Standards: IEC 60601-1, IEC 60601-1-2, IEC 60601-2-52



5-function ACP including
Lock-Out and Auto-Regression control



X-Ray Cassette



Head Angle Indicator /
Trend. Angle Indicator





Bed-end Extension 20cm

ES-12DF

Electric Hospital Bed

Standards: IEC 60601-1, IEC 60601-1-2, IEC 60601-2-52



Chair Position



Bed Height Indicator



Safety Standard
IEC 60601-2-52





Bed-end Extension 20cm

ES-11DF

Electric Hospital Bed

Standards: IEC 60601-1, IEC 60601-1-2, IEC 60601-2-52



CPR Quick Release (Option)

Attendant Control Panel (ACP)



Trend. / Rev. Trend. $\pm 12^\circ$ (Option)



CPR Quick Release (Option)



Drainage Bag Hanger



Bed-end Extension 10cm



Mattress Stopper



Specification	ES-12DW	ES-12HD	ES-12DF	ES-11DF
Overall bed length	L : 2220mm	L : 2220mm	L : 2220mm	L : 2120mm
Overall bed width (Siderails down/up)	W : 950mm 1000mm	W : 950mm 1000mm	W : 950mm 1000mm	W : 950mm 1000mm
Overall mattress base length	L : 2020mm	L : 2020mm	L : 2020mm	L : 1920mm
Overall mattress base width	W : 900mm	W : 900mm	W : 900mm	W : 900mm
Overall bed height	H : 455~755mm	H : 455~755mm	H : 455~755mm	H : 455~755mm
Head Elevation	0~75°	0~75°	0~75°	0~75°
Knee Elevation	0~30°	0~30°	0~30°	0~30°
Trend./Rev. Trend.	±12°	±12°	Optional	Optional
Work Load	230kg	230kg	230kg	230kg
Castor Diameter	5"	5"	5"	5"
Actuator	4	4	3/4	3
Volts	AC 100-240V	AC 100-240V	AC 100-240V	AC 100-240V
Frequency	50/60Hz	50/60Hz	50/60Hz	50/60Hz
Backup Battery	DC 24V 2.0Ah	DC 24V 2.0Ah	DC 24V 2.0Ah	DC 24V 2.0Ah
Actuator/Control Box	IP66	IP66	IP66	IP66
Safety Standard	IEC 60601-2-52	IEC 60601-2-52	IEC 60601-2-52	IEC 60601-2-52

Applicable mattress size (mm)	L2000xW880xH120	L2000xW880xH120	L2000xW880xH120	L1900xW880xH120

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree


Features	ES-12DW	ES-12HD	ES-12DF	ES-11DF
X-Ray Cassette	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Outer/Inner Control Panel	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Head Angle Indicator	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Trend./Rev. Trend. Angle Indicator	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
CPR Quick Release	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Roller Bumper	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Attendant Control Panel	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Remote Controller	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Weight Scale	<input checked="" type="checkbox"/>			
Lock-Out Function	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Auto Regression	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Backup Manual Control Module			<input type="checkbox"/>	<input type="checkbox"/>
Bed-end Extension	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Chair Position	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Trend./Rev. Trend.	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Dust-Proof Central Braking/Steering Module	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
IV Pole Hole	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Drainage Bags Hanger	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Mattress Stopper	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Height Caution Indicator	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Night Light	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
LCD Screen Panel Remote Controller	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Standard Option


09.96 series




ES-09DW



Head Elevation
0~85°



Knee Elevation
0~22°



Bed Elevation



Auto-chair



Trend. $\pm 12^\circ$



Rev. Trend. $\pm 12^\circ$



Head Elevation
indicator




Safety Standard



225KG

Work Load



4PCS

Actuator




Backup Battery




IP66

Actuator/Control
Box


ES-96HD



Head Elevation
0~85°




Knee Elevation
0~32°



Bed Elevation



Auto-chair



Trend. $\pm 12^\circ$



Rev. Trend. $\pm 12^\circ$



Trend.
Angle Indicator



Head
Angle Indicator




Safety Standard



225KG

Work Load



4PCS

Actuator



Backup Battery





IP66

Actuator/Control
Box

ES-09DW

Electric Hospital Bed
with Weight Scale




Standards: IEC 60601-1-2



ES-96HD

Electric Hospital Bed

Standards: IEC 60601-1-2, IEC 60601-2-38



Attendant Control Panel (ACP)



Drop Down Siderail



Mattress Stopper



Specification	ES-09DW	ES-96HD
Overall bed length	L : 2250mm	L : 2250mm
Overall bed width	W : 940mm	W : 940mm
(Siderails down/up)	1060mm	1060mm
Overall mattress base length	L : 2020mm	L : 2020mm
Overall mattress base width	W : 925mm	W : 925mm
Overall bed height	H : 455~755mm	H : 455~755mm
Head Elevation	0~85°	0~85°
Knee Elevation	0~22°	0~32°
Trend. / Rev. Trend.	±12°	±12°
Work Load	225kg	225kg
Castor Diameter	5"	5"
Actuator	4	4
Volts & Frequency	AC 100-240V ; 50/60Hz	AC 100-240V ; 50/60Hz
Backup Battery	DC 24V 2.0Ah	DC 24V 2.0Ah
Actuator/Control Box	IP66	IP66
Safety Standard	IEC 60601-2-38	IEC 60601-2-38
Applicable mattress size (mm)	L1980×W890×H120	L1980×W890×H100

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree


Features	ES-09DW	ES-96HD
X-Ray Cassette	<input checked="" type="checkbox"/>	
Outer/Inner Control Panel	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Backrest, Trend./Rev. Trend. Angle Indicator	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
CPR Quick Release	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Bumper	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Attendant Control Panel	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Electric CPR	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Remote Controller	<input type="checkbox"/>	<input type="checkbox"/>
Weight Scale	<input checked="" type="checkbox"/>	
Lock-Out Function	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Auto Regression	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Lockable Castor		<input checked="" type="checkbox"/>
Dust-Proof Central Braking/Steering Module	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Dust-Proof Lockable Caster		<input type="checkbox"/>
IV Pole Hole	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Drainage Bag Hanger	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Mattress Stopper	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

Standard Option

05.99 series




ES-05HD / ES-05HDS



Head Elevation
0 ~ 85°



Knee Elevation
0 ~ 32°



Bed Elevation



Auto-chair




Safety Standard



225KG

Work Load



3PCS

Actuator




Backup Battery



IP66

Actuator/Control
Box


MS-05HS



Head Elevation
0 ~ 85°



Knee Elevation
0 ~ 32°




Bed Elevation



225KG

Work Load

ES-99HDS



Head Elevation
0 ~ 85°



Knee Elevation
0 ~ 32°




Bed Elevation




225KG

Work Load



3PCS

Actuator



Backup Battery



Auto-chair



Safety Standard





IP66

Actuator/Control
Box

ES-05HD

Electric Hospital Bed

Standards: IEC 60601-1-2, IEC 60601-2-38



Attendant Control Panel (ACP)



Swivel Siderail




Backrest Angle Indicator (Option)

ES-05HDS

Electric Hospital Bed

Standards: IEC 60601-1-2, IEC 60601-2-38



Hospital Bed



Siderails up Bed width 100cm fit to small ward.



Detachable Head/Foot Board



IV Pole Hole





Drainage Bag Hanger



Mattress Stopper

MS-05HS



Manual Hospital Bed




ES-99HDS

Electric Hospital Bed

Standard: IEC 60601-1-2, IEC 60601-2-38



Manual Crank Screw drive



Mattress Stopper



Swivel Siderail



IV Pole Hole



Restraint Belt Hole



Specification	ES-05HD	ES-05HDS	MS-05HS	ES-99HDS
Overall bed length	L : 2170mm	L : 2115mm	L : 2170mm	L : 2115mm
Overall bed width (Siderails down/up)	W : 940mm 1090mm	W : 910mm 1000mm	W : 940mm 1090mm	W : 930mm 1090mm
Overall mattress base length	L : 1940mm	L : 1940mm	L : 1940mm	L : 1940mm
Overall mattress base width	W : 925mm	W : 845mm	W : 925mm	W : 925mm
Overall bed height (Central Braking/Normal Castor)	H : 455~755mm H : 395~695mm	H : 455~755mm H : 395~695mm	H : 455~755mm H : 395~695mm	H : 455~755mm H : 395~695mm
Head Elevation	0~85°	0~85°	0~85°	0~85°
Knee Elevation	0~32°	0~32°	0~32°	0~32°
Trend./Rev. Trend.	None	None	None	None
Work Load	225kg	225kg	225kg	225kg
Caster Diameter	5"	5"	5"	5"
Actuator	3	3	Manual Crank	3
Volts & Frequency	AC 100-240V ; 50/60Hz	AC 100-240V ; 50/60Hz	None	AC 100-240V ; 50/60Hz
Backup Battery	DC 24V 2.0Ah	DC 24V 2.0Ah	None	DC 24V 2.0Ah
Actuator/Control Box	IP66	IP66	None	IP66
Safety Standard	IEC 60601-2-38	IEC 60601-2-38	None	IEC 60601-2-38
Applicable mattress size (mm)	L1900xW890xH100	L1930xW840xH100	L1900xW890xH100	L1900xW890xH100

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree

Features	ES-05HD	ES-05HDS	MS-05HS	ES-99HDS
Outer/Inner Siderail Control	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Backrest, Trend./Rev. Trend. Angle Indicator	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
CPR Quick Release	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Bumper	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Attendant Control Panel	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Electric CPR	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Remote Controller	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Lock-Out Function	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Auto Regression	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pneumatic Cylinder assisted Hi-Lo Mechanism	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Manual Trend./Rev. Trend.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Dust-Proof Central Braking/Steering Module	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Dust-Proof Lockable Castor	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Lockable Castor	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
IV Pole Hole	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Drainage Bag Hanger	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Mattress Stopper	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

Standard Option

02.16 series




ES-02JH



ES-16N





MS-16NS



ES-02JH

Electric Nursing Bed

Standards: IEC 60601-1-2, IEC 60601-2-52



Rotary Armrest Siderail (Option)



Locking mechanism for fixed angles.



Detachable Head/Foot Board




Siderail Holder



Remote Controller



IV Pole Hole




Mattress Stopper

ES-16N series

Electric Nursing Bed

Standards: IEC 60601-1-2, IEC 60601-2-52

3/4 Length Siderail



Hospital Bed



3/4 Length Siderail + P-Type Siderail



Attendant Control Panel (Option)



Remote Controller



Detachable Head/Foot Board




Mattress Stopper

MS-16NS

Manual Hospital Bed

Standards: IEC 60601-1-2, IEC 60601-2-52



3/4 Length Siderail + P-Type Siderail



Manual Screw drive



IV Pole Hole



Detachable Head/Foot Board



Mattress Stopper



Specification	ES-02JH	ES-16N	MS-16NS
Overall bed length	L : 2130mm	L : 2115mm	L : 2115mm
Overall bed width	W : 1000mm	W : 1000mm	W : 1000mm
Overall mattress base length	L : 1910mm	L : 1910mm	L : 1910mm
Overall mattress base width	W : 845mm	W : 845mm	W : 845mm
Overall bed height	H : 320~630mm	H : 350~660mm	H : 350~660mm
Head Elevation	0~75°	0~85°	0~85°
Knee Elevation	0~30°	0~40°	0~40°
Work Load	170kg	225kg	225kg
Castor Diameter	4"	5"	5"
Actuator	3	3	Manual Crank
Volts	AC 120V	AC 120V	None
Frequency	50/60Hz	50/60Hz	None
Backup Battery	Optional	Optional	None
Actuator/Control Box	IP54	IP54	None
Safety Standard	IEC 60601-2-52	IEC 60601-2-52	IEC 60601-2-52

Applicable mattress size (mm)	L1930xW840xH100	L1930xW840xH100	L1930xW840xH100

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree (ES-02JH ±3°)


Features	ES-02JH	ES-16N	MS-16NS
Remote Controller	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Siderail	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Bumper	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Attendant Control Panel	<input type="checkbox"/>	<input type="checkbox"/>	
Lock-Out Function		<input type="checkbox"/>	
Auto Regression		<input type="checkbox"/>	
Pneumatic Cylinder assisted Hi-Lo Mechanism			<input checked="" type="checkbox"/>
Dust-Proof Central Braking /Steering Module	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Dust-Proof Lockable Castor	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Lockable Castor		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
IV Pole Hole	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Drainage Bag Hanger	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Mattress Stopper	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

Standard Option


Stretcher series




JE-350 & JE-300



Head Elevation
0 ~ 90°



Bed Elevation



Trend. $\pm 16^\circ$



Rev. Trend. $\pm 16^\circ$



225KG

Work Load

JE-200 & JE-150



Head Elevation
0 ~ 90°



Knee Elevation
0 ~ 30°



Bed Elevation



Trend. $\pm 18^\circ$



Rev. Trend. $\pm 18^\circ$






225KG

Work Load

JE-350 Hydraulic Stretcher



The Fifth Wheel reduces the floor surface required when turning the bed.



unit: mm



X-Ray Cassette



Whole Body X-Ray Cassette System



Functional Pedal






Monitor Tray

JE-300 Hydraulic Stretcher



The Fifth Wheel reduces the floor surface required when turning the bed.



unit: mm



Back Gas Spring



Siderail Lock Handgrip



Oxygen Cylinder



3-section IV Pole




Foot Gas Spring
(Option)

JE-200

Hydraulic Stretcher



The Fifth Wheel reduces the floor surface required when turning the bed.



X-Ray Cassette (Option)



Retractable Push Handgrip




Oxygen Cylinder



Paper Roll Holder

JE-150

Hydraulic Stretcher



Lifting Rod



Knee Elevation Crank



Collapsible Siderail



Siderail Lock Handgrip



Dust-Proof Central Braking Steering Module



Specification	JE-350	JE-300	JE-200	JE-150	Features	JE-350	JE-300	JE-200	JE-150
Overall bed length	L : 2120mm	L : 2120mm	L : 2100mm	L : 2100mm	X-Ray	✓		□	
Overall bed width	W : 790mm	W : 790mm	W : 825mm	W : 825mm	Full-Length X-Ray Capacity	✓			
(Siderails down/up)	W : 830mm	W : 830mm	W : 840mm	W : 840mm	Outward Siderail	✓	✓		
Overall mattress base length	L : 1850mm	L : 1850mm	L : 1900mm	L : 1900mm	Collapsible Siderail			✓	✓
Overall mattress base width	W : 635mm	W : 635mm	W : 740mm	W : 740mm	Knee Elevation Handgrip	✓	✓	✓	✓
Overall bed height	H : 590~890mm	H : 590~890mm	H : 585~935mm	H : 585~935mm	Roller Bumper	✓	✓		
Head Elevation	0~90°	0~90°	0~90°	0~90°	Bumper			✓	✓
Knee Elevation	None	None	0~22°	0~22°	Dust-Proof Central Braking/Steering Module	✓	✓	✓	✓
Trend./Rev. Trend.	16°	16°	18°	18°	The Fifth Wheel Module	✓	✓	✓	
Work Load	225kg	225kg	225kg	225kg	Oxygen Bottle Holder	✓	✓	□	□
Frame Material	Power Coated Mild Steel	Power Coated Mild Steel	Power Coated Mild Steel	Power Coated Mild Steel	Storage Tray	✓	✓	✓	✓
Castor Diameter	8"	8"	8"	8"	IV Pole	✓	✓	✓	✓
Transmission	Backrest: Pneumatic Cylinder Knee Joint: Hand-Held Bed Hi-Lo: Hydraulic Cylinder	Backrest: Pneumatic Cylinder Knee Joint: Hand-Held Bed Hi-Lo: Hydraulic Cylinder	Backrest: Pneumatic Cylinder Knee Joint: Hand-Held; Crank Bed Hi-Lo: Hydraulic Cylinder	Backrest: Pneumatic Cylinder Knee Joint: Hand-Held; Crank Bed Hi-Lo: Hydraulic Cylinder	Drainage Bag Hanger	✓	✓	□	□
Applicable mattress size(mm)	L1820xW620xH75	L1820xW620xH75	L1850xW720xH100	L1850xW720xH100	Restraint Belt Hole	□	□	□	□

Standard
 Option

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree

JE-130

3 Manual Cranks
Emergency Stretcher



JE-120

2 Manual Cranks
Emergency Stretcher



JE-100

Emergency Stretcher



Specification	JE-130	JE-120	JE-100
Overall bed length	L : 1900mm	L : 1900mm	L : 1900mm
Overall bed width	W : 750mm	W : 750mm	W : 750mm
Overall mattress base length	L : 1810mm	L : 1810mm	L : 1810mm
Overall mattress base width	W : 550mm	W : 550mm	W : 550mm
Overall bed height	H : 550~750mm	H : 650mm	H : 650mm
Head Elevation	0~80°	0~80°	0~80°
Knee Elevation	0~35°	0~35°	None
Work Load	150kg	150kg	150kg
Frame Material	SUS304	SUS304	SUS304
Castor Diameter	6"	6"	6"
Transmission	Crank	Crank	Hand-Hold
Mattress	2" High density PU foam	2" High density PU foam	2" High density PU foam

Features	JE-130	JE-120	JE-100
Collapsible Siderail	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Lockable Castor	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Oxygen Bottle Holder	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Storage Tray	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
IV Pole	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

Standard


Applicable mattress size (mm)	L1800xW550xH50	L1800xW550xH50	L1800xW550xH50

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree

Pediatric Bed series




JC-012



Head Elevation
0~80°



Knee Elevation
0~40°




Bed Elevation




170KG

Work Load



3PCS

Actuator




Backup Battery



IP66

Actuator/Control
Box

JC-011



Head Elevation
0~80°



Knee Elevation
0~40°



90KG

Work Load

JC-010



Head Elevation
0~80°



90KG

Work Load


JC-012

Electric Pediatric Bed



Safety Siderail Handle

Attendant Control Panel (ACP)



unit: mm



Attendant Control Panel (ACP)



Safety Siderail Handle




Remote Controller



Dust-Proof Lockable Castor

JC-011

Manual Pediatric Bed



Manual 2 Cranks




Solid Structure



Safety Siderail Handle

JC-010

Manual Pediatric Bed



JC-031 Infant Crib

Size : L700 x W440 x H910mm
 ±20mm tolerance in all dimensions
 Frame : ABS Plastic, Acrylic Basin
 Basin Angle: 0~20°



JC-030 Infant Crib

Size : L680 x W380 x H760mm
 ±20mm tolerance in all dimensions
 Frame : SUS304, Acrylic Baby Basin



Specification	JC-012	JC-011	JC-010
Overall bed length	L : 2050mm	L : 2000mm	L : 1600mm
Overall bed width	W : 1000mm	W : 816mm	W : 816mm
Overall mattress base length	L : 1870mm	L : 1944mm	L : 1545mm
Overall mattress base width	W : 845mm	W : 756mm	W : 758mm
Overall bed height	H : 650~850mm	H : 763mm	H : 764mm
Head Elevation	0~80°	0~80°	0~80°
Knee Elevation	0~40°	0~40°	None
Work Load	170kg	90kg	90kg
Frame Material	SUS304+Powder Coated	SUS304+Powder Coated	SUS304+Powder Coated
Castor Diameter	5"	4"	4"
Actuator	3	Manual (2 Cranks)	Manual (1 Crank)
Volts & Frequency	AC 100-240V ; 50/60Hz	None	None
Backup Battery	DC 24V 2.0Ah	None	None
Actuator/Control Box	IP66	None	None

Applicable mattress size (mm)	JC-012	JC-011	JC-010
	L1850xW830xH30	L1920xW720xH30	L1440xW780xH76

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree

Features	JC-012	JC-011	JC-010
Remote Controller	<input checked="" type="checkbox"/>		
Bumper	<input checked="" type="checkbox"/>		
Attendant Control Panel	<input checked="" type="checkbox"/>		
Electric CPR	<input checked="" type="checkbox"/>		
Lock-Out Function	<input checked="" type="checkbox"/>		
Auto Regression	<input checked="" type="checkbox"/>		
Dust-Proof Central Braking/Steering Module	<input checked="" type="checkbox"/>		
Dust-Proof Lockable Castor	<input type="checkbox"/>		
Lockable Castor		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
IV Pole Hole	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

Standard Option

18 series



ES-18UDS / ES-17UDS



Head Elevation
0 ~ 70°



Knee Elevation
0 ~ 30°




Bed Elevation



Auto-chair



170KG
Work Load



3PCS
Actuator



Backup
Battery Box
(Option)




IP54
Actuator/Control
Box




Shear Stress
Bi-directional Release

ES-18UDS

Ultra Low Bed



LO: 23.5cm



Shear Stress
Bi-directional Release



Dust-Proof Central Braking
Steering Module




Detachable
Head/Foot Board



Mattress Stopper

ES-17UDS

Ultra Low Bed



Home Care Bed



Shear Stress Bi-directional Release



Detachable Head/Foot Board



Dust-Proof Central Braking Steering Module



3/4 Length Siderail, Rotary Armrest Siderail, Plug-in Siderail

18 / 17 series

Ultra Low Design

Ultra low at 23.5 cm



ES-18UDS



ES-17UDS

Shear Stress Bi-directional Release

The hip and back boards decline to minimize shear stress and ease pressure on the back and abdomen.



More space between back board and hip board.





Specification	ES-18UDS	ES-17UDS	Features	ES-18UDS	ES-17UDS
Overall bed length	L : 2080mm	L : 2100mm	Remote Controller	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Overall bed width	W : 1050mm	W : 950mm	Rotary Siderail	<input type="checkbox"/>	<input type="checkbox"/>
Overall mattress base length	L : 1915mm	L : 1915mm	3/4 Length Siderail	<input type="checkbox"/>	<input type="checkbox"/>
Overall mattress base width	W : 910mm	W : 830mm	Plug-in Siderail	<input type="checkbox"/>	<input type="checkbox"/>
Overall bed height	4"Lockable Caster 235~535mm	4"Lockable Caster 235~535mm	Detachable Head/Foot Board	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Head Elevation	0~70°	0~70°	Dust-Proof Central Braking / Steering Module	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Knee Elevation	0~30°	0~30°	Dust-Proof Lockable Castor	<input type="checkbox"/>	<input type="checkbox"/>
Work Load	170kg	170kg	Lockable Castor	<input type="checkbox"/>	<input type="checkbox"/>
Castor Diameter	4"	4"	IV Pole Hole	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Actuator	3	3	Mattress Stopper	<input checked="" type="checkbox"/>	
Volts & Frequency	AC110V ; 50/60Hz	AC110V ; 50/60Hz	Backup Battery (option) 9V Alkaline Battery x3	<input type="checkbox"/>	<input type="checkbox"/>
Actuator/Control Box	IP54	IP54			

Standard Option


Applicable mattress size (mm)	L1900×W890×H100	L1830×W830×H100
-------------------------------	-----------------	-----------------

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree


EN series



EN-3M





EN-2M



EN-3M

Foldable Electric Home Care Bed



Trend./Rev. Trend. $\pm 11^\circ$



Chair Position



Wooden Head/Foot Board





3/4 Length Siderail, Rotary Armrest Siderail, Plug-in Siderail

EN-2M

Assembly Electric Home Care Bed

5 Cartons Pack



Trend./Rev. Trend. $\pm 11^\circ$




Chair Position



Wooden Head/Foot Board



Dust-Proof Lockable Castor



Mattress Stopper


3 Siderails



Rotary Armrest Siderail
Siderail installs on bedside, caregivers could adjust angles upon patient's need.



Plug-in Siderail
Siderail installs on bedside or foot-end.



3/4 Length Siderail

Siderails for below models

ES-02JH | EN-3M | EN-2M | ES-18UDS | ES-17UDS



ES-02JH



EN-3M



EN-2M



ES-18UDS



Specification	EN-3M	EN-2M	Features	EN-3M	EN-2M
Overall bed length	L : 2140mm	L : 2140mm	Chair Position	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Overall bed width	W : 970mm	W : 1070mm	Foldable Type	<input checked="" type="checkbox"/>	
Overall mattress base length	L : 1935mm	L : 1900mm	Assembly Type		<input checked="" type="checkbox"/>
Overall mattress base width	W : 775mm	W : 910mm	Remote Controller	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Overall bed height	4"Lockable Caster 300~700mm	4"Lockable Caster 300~700mm	Rotary Armrest Siderail	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Head Elevation	0~70°	0~75°	3/4 Length Siderail	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Knee Elevation	0~28°	0~15°	Plug-in Siderail	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Trend./Rev. Trend.	11°	11°	Dust-Proof Lockable Castor	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Work Load	170kg	173.5kg	Lockable Castor	<input checked="" type="checkbox"/>	
Castor Diameter	4"	4"	IV Pole Hole	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Actuator	4	3	Mattress Stopper		<input checked="" type="checkbox"/>
Volts & Frequency	AC110 ; 50/60Hz	AC110V ; 50/60Hz	Backup Battery (option) 9V Alkaline Battery x3	<input type="checkbox"/>	<input type="checkbox"/>
Actuator/Control Box	IP54	IP54			

Standard Option


Applicable mattress size (mm)	L1980xW820xH150	L1980xW890xH100
-------------------------------	-----------------	-----------------

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree

08 series




ES-08FDS



Head Elevation
0~80°



Knee Elevation
0~40°



Bed Elevation




Auto-chair



170KG

Work Load



3PCS


Actuator



IP42

Actuator/Control
Box

ES-08FSS



Head Elevation
0~85°



Knee Elevation
0~40°



Bed Elevation




Auto-chair



170KG

Work Load



3PCS


Actuator



IP42

Actuator/Control
Box


MS-08FS



Head Elevation
0~80°



Knee Elevation
0~40°



Bed Elevation




170KG

Work Load


ES-08FDS

Electric Hospital Bed

3/4 Length Siderail



4" castor
H: 420~620
5" castor
H: 430~630



unit: mm




3/4 Length Siderail



Remote Controller




Drainage Bag Hanger



Plug-in
Head/Foot Board



Dust-Proof Central Braking
Steering Module



Mattress Stopper


ES-08FSS

Electric Hospital Bed
(Metal Frame)





Remote Controller

3/4 Length Siderail



4" castor
H: 420~620
5" castor
H: 430~630



unit: mm



Metal Frame




Remote Controller



IV Pole Hole




Lockable Castor



Bed-end Extension
(Option)

MS-08FS / MF-08FS


3 Manual Cranks / 2 Manual Cranks



4" castor
H: 420~620

5" castor
H: 430~630

MF-08FS H: 450



unit: mm

3/4 Length Siderail




Manual Crank



Drainage Bag Hanger



3/4 Length Siderail



Mattress Stopper



Specification	ES-08FDS	ES-08FSS	MS-08FS	MF-08FS	Features	ES-08FDS	ES-08FSS	MS-08FS	MF-08FS
Overall bed length	L : 2010mm	L : 1960mm	L : 2010mm	L : 2010mm	Remote Controller	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
Overall bed width	W : 940mm	W : 930mm	W : 940mm	W : 940mm	3/4 Length Siderail	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Overall mattress base length	L : 1860mm	L : 1770mm	L : 1860mm	L : 1860mm	Dust-Proof Central Braking/Steering Module	<input type="checkbox"/>	<input type="checkbox"/>		
Overall mattress base width	W : 845mm	W : 815mm	W : 845mm	W : 845mm	Dust-Proof Lockable Castor	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Overall bed height	H:420~620mm (4" castor) H:430~630mm (5" castor)	H:420~620mm (4" castor) H:430~630mm (5" castor)	H:420~620mm (4" castor) H:430~630mm (5" castor)	H : 450mm	Lockable Castor	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Head Elevation	0~80°	0~85°	0~80°	0~80°	Bed-Head Extension		<input checked="" type="checkbox"/>		
Knee Elevation	0~40°	0~40°	0~40°	0~40°	IV Pole Hole	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Work Load	170kg	170kg	170kg	170kg	Drainage Bag Hanger	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Castor Diameter	5" or 4"	5" or 4"	5" or 4"	5"	Mattress Stopper	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Actuator	3	3	Manual Crank	Manual Crank	<input checked="" type="checkbox"/> Standard <input type="checkbox"/> Option				
Volts	AC 115V/220V	AC 115V/220V	none	none					
Frequency	50/60Hz	50/60Hz	none	none					
Backup Battery	none	none	none	none					
Actuator/Control Box	IP42	IP42	none	none					
Safety Standard	none	none	none	none					

Applicable mattress size (mm)	ES-08FDS	ES-08FSS	MS-08FS	MF-08FS
	L1850xW840xH100	L1800xW810xH100	L1850xW840xH100	L1850xW840xH100

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree


Mattress

JS-301, JS-302, JS-303

Mattress cover: Nylon



■ BS



Ventilation



Waterproof




Flame Retardant



Anti Bacterial

■ CA



Ventilation




Waterproof




Flame Retardant


High-Density Foam



JS-301



JS-302



JS-303

Specification	JS-301		JS-302		JS-303	
Model	JS-301BS	JS-301CA	JS-302BS	JS-302CA	JS-303BS	JS-303CA
Mattress cover	■	■	■	■	■	■
Fire safety standard	BS-7177	CA-117	BS-7177	CA-117	BS-7177	CA-117
Anti-bacterial	✓		✓		✓	

JT-100

Orthopedic Traction Frame



Triangle Ring Set



Pulley Mounting Set



Foot Support Unit



Support Rod Set



1. The traction frame is suitable for all series of Joson-Care's bed and installed in four IV Pole holder.
2. The bottom of the frame is made by solid stainless steel bars.
3. Switching the wrench could quickly adjust the frame, it's tools free to set the frame.
4. The traction frame includes one set of triangular ring and pulley.
5. The octagonal pipes are made of aluminum.
6. Two sets of reinforcing units are helping to reinforce the traction frame.
7. Adjustable height and length of triangle ring set can avoid patient being impacted while position changed.

IV Stand



JR-016 IV Stand

Size : H1630~2300mm
Frame : Stainless 304#



JR-018 Stackable IV Stand

Size : 1335~2169mm
Frame : Chrome Plated steel



K5000-MS 2 Claws IV Pole

Size : 855~1510mm
Frame : Chrome Plated steel



K4000-MS 4 Claws IV Pole

Size : 1020~1760mm
Frame : Chrome Plated steel

JA-107

Emergency Cart

Overall Size : L735 x W560 x H990mm
Table Surface Size : 380 x 817 mm
Drawers size : 60 x 9 x 15 : ±1 cm




Laptop Drawer
(390 x 380mm)



Trash Bin(9 liters)
(W277 x H182 x D390mm)



Cardiac Board
(L700 x W400 x D15mm)



Disposable Lock

- Frame: Steel
- Cover: ABS
- A buckle fixed design which prevent drawers sliding out.
- Key and disposable lock available.
- Vertical and horizontal divider included.
- Flexible confuguration design.
- 2 sets divider quipe in both 3" and 6" drawers.

Accessories Option



Storage Box
(W600 x H304 x D124mm)



IV Pole
(1020~1760mm)



Suction Tube Holder
(W90 x H480 x D90mm)



Oxygen Tank Holder
(W125 x H370 x D175mm)



Pedal Trashbin(30 liters)
(W311 x H406 x D500mm)



Hand Sanitizer Holder
(Including hand sanitizer bottle)
(W75 x H105 x D77mm)



Syringe Collection Box Holder
(W130 x H150 x D90mm)



Plastic Bin Holder
(W480 x H120 x D120mm)



Smesh Bin Holder
(W488 x H125 x D125mm)



Glove Dispenser

UD Medical Cart



Medical Cart



JO-103 Mayo Table

Size : L485 x W325 x H800~1200mm
Frame : Stainless 304#



JA-108 Dressing Trolley

Size : L600 x W450 x H850mm
Frame : Stainless 304#

JA-111 Utility Cart

Size : L600 x W450 x H850mm
Frame : Stainless 304#



JF-014

Bedside Cabinet



Size : L460 x W460 x H860mm
 Color : Blue / Pink
 Material: ABS

Outpatient Equipment



JR-017 Screen

Size : W750 x H1800mm
 Frame : Stainless 304#
 Cloth : 65% tricot & 35% rayon

JO-104 Pedal

Single Size : L400 x W300 x H250mm
 Double Size : L400 x W600 x H440mm
 Frame : Stainless 304#
 Surface : Anti-slip Plastic Board



Overbed Table

JR-032 Overbed Table

Size : L815 x W403 x H820~1225mm
 Frame : Powder coated steel
 Table Surface : flame-retardant laminate
 Max. Loading : 20KG



Gas spring driven design



2" lockable castor*4
 2 with brake



JR-011 Overbed Table

Size : L765 x W385 x H720~1120mm
 Frame : Powder coated steel
 Table Surface : flame-retardant laminate

JR-051 Overbed Table

Size : L845 x W400 x H780~1050mm
 Package : L880 x W465 x H165mm
 Frame : Powder coated steel
 Table Surface : flame-retardant laminate
 Max. Loading : 20KG



Y0307-SY Overbed Table

Size : L890~1060 x W320 x H30mm
 Frame : ABS plastic
 Tube : Aluminum
 Hooks on the backside.

Doctor Chair



JO-107 Doctor Chair

Size : Ø380mm · H510~660mm
 Seat: PU Foam
 Frame: Aluminum



JO-109 Doctor Chair

Size : Ø300mm · H500~670mm
 Frame : Stainless 304#



JO-108 Doctor Chair

Size : Ø380mm · H510~660mm
 Seat: PU Foam
 Frame: Aluminum

Caregiver Recliner



JD-010 Caregiver Recliner

Size : 600 x 800 x 950mm
 600 x 1850mm (Extended)

Frame : Stainless 304#
 Cushion : High density foam
 Cover : Synthetic Leather





Joson-Care

We serve all your care



Joson-Care Enterprise Co., Ltd.

No.11, Aly. 77, Ln. 180, Hushan St., Guishan Dist., Taoyuan City 333024, Taiwan

☎ +886-3-329-0925

✉ overseas@joson-care.com.tw

🌐 www.joson-care.com

Taichung Office

1F., No. 280, Chengdu Rd., Xitun Dist.,
Taichung City 407403, Taiwan

TEL : +886-4-2316-0743

Tainan Office

No. 732, Sec. 1, Zhongzheng Rd., Rende Dist.,
Tainan City 717010, Taiwan

TEL : +886-6-266-2809



Document No.: JCQB-031
Release Date: 2024.06.01
Version: L