



# Joson-Care 2024




Safety Quality Sustainability






Business  
Concept




History




Established in 2000, Joson-Care Co., Ltd. is Taiwan's premier manufacturer of hospital beds and ward equipment. We uphold rigorous quality control standards to deliver exceptional accuracy, durability, and reliability, while maximizing the use of renewable resources.





So far, we have sold our products to more than 40 countries. Such as U.S.A, Japan, France, South Korea, India, Iran, Panama, Thailand, Turkey and U.A.E.


Joson-Care's manufacturing system is meticulously designed to adhere to ISO 13485 and Quality Management System (QMS) standards. This rigorous approach ensures that our processes consistently deliver high-quality products, precisely meeting and often exceeding customer expectations for safety, performance, and reliability.




# Stretcher series




## JE-350 & JE-300



Head Elevation  
0 ~ 90°



Bed Elevation



Trend.  $\pm 16^\circ$



Rev. Trend.  $\pm 16^\circ$



225KG

Work Load

## JE-200 & JE-150



Head Elevation  
0 ~ 90°



Knee Elevation  
0 ~ 30°



Bed Elevation



Trend.  $\pm 18^\circ$



Rev. Trend.  $\pm 18^\circ$



225KG





Work Load

# JE-350

## Hydraulic Stretcher



The Fifth Wheel reduces the floor surface required when turning the bed.



unit: mm



X-Ray Cassette



Whole Body X-Ray Cassette System



Functional Pedal





Monitor Tray

# JE-300 Hydraulic Stretcher



The Fifth Wheel reduces the floor surface required when turning the bed.



Back Gas Spring



Siderail Lock Handgrip



Oxygen Cylinder



3-section IV Pole




Foot Gas Spring  
(Option)

# JE-200

## Hydraulic Stretcher



The Fifth Wheel reduces the floor surface required when turning the bed.



X-Ray Cassette (Option)



Retractable Push Handgrip



Oxygen Cylinder




Paper Roll Holder



# JE-150

## Hydraulic Stretcher



Lifting Rod



Knee Elevation Crank



Collapsible Siderail



Siderail Lock Handgrip



Dust-Proof Central Braking Steering Module



Specification	JE-350	JE-300	JE-200	JE-150	Features	JE-350	JE-300	JE-200	JE-150
Overall bed length	L : 2120mm	L : 2120mm	L : 2100mm	L : 2100mm	X-Ray	<input checked="" type="checkbox"/>		<input type="checkbox"/>	
Overall bed width	W : 790mm	W : 790mm	W : 825mm	W : 825mm	Full-Length X-Ray Capacity	<input checked="" type="checkbox"/>			
(Siderails down/up)	W : 830mm	W : 830mm	W : 840mm	W : 840mm	Outward Siderail	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
Overall mattress base length	L : 1850mm	L : 1850mm	L : 1900mm	L : 1900mm	Collapsible Siderail			<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Overall mattress base width	W : 635mm	W : 635mm	W : 740mm	W : 740mm	Knee Elevation Handgrip	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Overall bed height	H : 590~890mm	H : 590~890mm	H : 585~935mm	H : 585~935mm	Roller Bumper	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
Head Elevation	0~90°	0~90°	0~90°	0~90°	Bumper			<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Knee Elevation	None	None	0~22°	0~22°	Dust-Proof Central Braking/Steering Module	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Trend./Rev. Trend.	16°	16°	18°	18°	The Fifth Wheel Module	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Work Load	225kg	225kg	225kg	225kg	Oxygen Bottle Holder	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Frame Material	Power Coated Mild Steel	Power Coated Mild Steel	Power Coated Mild Steel	Power Coated Mild Steel	Storage Tray	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Castor Diameter	8"	8"	8"	8"	IV Pole	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Transmission	Backrest: Pneumatic Cylinder Knee Joint: Hand-Held Bed Hi-Lo: Hydraulic Cylinder	Backrest: Pneumatic Cylinder Knee Joint: Hand-Held Bed Hi-Lo: Hydraulic Cylinder	Backrest: Pneumatic Cylinder Knee Joint: Hand-Held; Crank Bed Hi-Lo: Hydraulic Cylinder	Backrest: Pneumatic Cylinder Knee Joint: Hand-Held; Crank Bed Hi-Lo: Hydraulic Cylinder	Drainage Bag Hanger	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Applicable mattress size(mm)	L1820xW620xH75	L1820xW620xH75	L1850xW720xH100	L1850xW720xH100	Restraint Belt Hole	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Standard     Option

±20mm tolerance in all dimensions    ±30mm in siderail width    +5° -0° tolerance in angle degree

## JE-130

3 Manual Cranks  
Emergency Stretcher



## JE-120

2 Manual Cranks  
Emergency Stretcher



## JE-100

Emergency Stretcher



Specification	JE-130	JE-120	JE-100
Overall bed length	L : 1900mm	L : 1900mm	L : 1900mm
Overall bed width	W : 750mm	W : 750mm	W : 750mm
Overall mattress base length	L : 1810mm	L : 1810mm	L : 1810mm
Overall mattress base width	W : 550mm	W : 550mm	W : 550mm
Overall bed height	H : 550~750mm	H : 650mm	H : 650mm
Head Elevation	0~80°	0~80°	0~80°
Knee Elevation	0~35°	0~35°	None
Work Load	150kg	150kg	150kg
Frame Material	SUS304	SUS304	SUS304
Castor Diameter	6"	6"	6"
Transmission	Crank	Crank	Hand-Hold
Mattress	2" High density PU foam	2" High density PU foam	2" High density PU foam

Applicable mattress size (mm)	L1800xW550xH50	L1800xW550xH50	L1800xW550xH50
-------------------------------	----------------	----------------	----------------

±20mm tolerance in all dimensions    ±30mm in siderail width    +5° -0° tolerance in angle degree

Features	JE-130	JE-120	JE-100
Collapsible Siderail	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Lockable Castor	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Oxygen Bottle Holder	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Storage Tray	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
IV Pole	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

Standard



### Joson-Care Enterprise Co., Ltd.

No.11, Aly. 77, Ln. 180, Hushan St., Guishan Dist., Taoyuan City 333024, Taiwan

☎ +886-3-329-0925

✉ overseas@joson-care.com.tw

🌐 www.joson-care.com

#### Taichung Office

1F., No. 280, Chengdu Rd., Xitun Dist.,  
Taichung City 407403, Taiwan

TEL : +886-4-2316-0743

#### Tainan Office

No. 732, Sec. 1, Zhongzheng Rd., Rende Dist.,  
Tainan City 717010, Taiwan

TEL : +886-6-266-2809



Document No.: JCQB-031  
Release Date: 2024.06.01  
Version: L