



# Joson-Care 2024




Safety Quality Sustainability






Business  
Concept




History




Established in 2000, Joson-Care Co., Ltd. is Taiwan's premier manufacturer of hospital beds and ward equipment. We uphold rigorous quality control standards to deliver exceptional accuracy, durability, and reliability, while maximizing the use of renewable resources.





So far, we have sold our products to more than 40 countries. Such as U.S.A, Japan, France, South Korea, India, Iran, Panama, Thailand, Turkey and U.A.E.


Joson-Care's manufacturing system is meticulously designed to adhere to ISO 13485 and Quality Management System (QMS) standards. This rigorous approach ensures that our processes consistently deliver high-quality products, precisely meeting and often exceeding customer expectations for safety, performance, and reliability.



# 18 series



## ES-18UDS / ES-17UDS



Head Elevation  
0 ~ 70°



Knee Elevation  
0 ~ 30°



Bed Elevation




Auto-chair



170KG

Work Load



3PCS

Actuator



Backup  
Battery Box  
(Option)



IP54



Actuator/Control  
Box




Shear Stress  
Bi-directional Release

# ES-18UDS

## Ultra Low Bed



LO: 23.5cm



Shear Stress  
Bi-directional Release



Dust-Proof Central Braking  
Steering Module




Detachable  
Head/Foot Board



Mattress Stopper

# ES-17UDS

## Ultra Low Bed



Home Care Bed



Shear Stress Bi-directional Release



Detachable Head/Foot Board



Dust-Proof Central Braking Steering Module




3/4 Length Siderail, Rotary Armrest Siderail, Plug-in Siderail

## 18 / 17 series

Ultra Low Design

Ultra low at 23.5 cm



ES-18UDS



ES-17UDS

## Shear Stress Bi-directional Release

The hip and back boards decline to minimize shear stress and ease pressure on the back and abdomen.



More space between back board and hip board.







Specification	ES-18UDS	ES-17UDS	Features	ES-18UDS	ES-17UDS
Overall bed length	L : 2080mm	L : 2100mm	Remote Controller	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Overall bed width	W : 1050mm	W : 950mm	Rotary Siderail	<input type="checkbox"/>	<input type="checkbox"/>
Overall mattress base length	L : 1915mm	L : 1915mm	3/4 Length Siderail	<input type="checkbox"/>	<input type="checkbox"/>
Overall mattress base width	W : 910mm	W : 830mm	Plug-in Siderail	<input type="checkbox"/>	<input type="checkbox"/>
Overall bed height	4"Lockable Caster 235~535mm	4"Lockable Caster 235~535mm	Detachable Head/Foot Board	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Head Elevation	0~70°	0~70°	Dust-Proof Central Braking / Steering Module	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Knee Elevation	0~30°	0~30°	Dust-Proof Lockable Castor	<input type="checkbox"/>	<input type="checkbox"/>
Work Load	170kg	170kg	Lockable Castor	<input type="checkbox"/>	<input type="checkbox"/>
Castor Diameter	4"	4"	IV Pole Hole	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Actuator	3	3	Mattress Stopper	<input checked="" type="checkbox"/>	
Volts & Frequency	AC110V ; 50/60Hz	AC110V ; 50/60Hz	Backup Battery (option) 9V Alkaline Battery x3	<input type="checkbox"/>	<input type="checkbox"/>
Actuator/Control Box	IP54	IP54			

Standard     Option

Applicable mattress size (mm)	L1900×W890×H100	L1830×W830×H100
-------------------------------	-----------------	-----------------

±20mm tolerance in all dimensions    ±30mm in siderail width    +5° -0° tolerance in angle degree



### Joson-Care Enterprise Co., Ltd.

No.11, Aly. 77, Ln. 180, Hushan St., Guishan Dist., Taoyuan City 333024, Taiwan

☎ +886-3-329-0925

✉ overseas@joson-care.com.tw

🌐 www.joson-care.com

#### Taichung Office

1F., No. 280, Chengdu Rd., Xitun Dist.,  
Taichung City 407403, Taiwan

TEL : +886-4-2316-0743

#### Tainan Office

No. 732, Sec. 1, Zhongzheng Rd., Rende Dist.,  
Tainan City 717010, Taiwan

TEL : +886-6-266-2809



Document No.: JCQB-031  
Release Date: 2024.06.01  
Version: L