



Joson-Care 2024




Safety Quality Sustainability






Business
Concept




History




Established in 2000, Joson-Care Co., Ltd. is Taiwan's premier manufacturer of hospital beds and ward equipment. We uphold rigorous quality control standards to deliver exceptional accuracy, durability, and reliability, while maximizing the use of renewable resources.





So far, we have sold our products to more than 40 countries. Such as U.S.A, Japan, France, South Korea, India, Iran, Panama, Thailand, Turkey and U.A.E.


Joson-Care's manufacturing system is meticulously designed to adhere to ISO 13485 and Quality Management System (QMS) standards. This rigorous approach ensures that our processes consistently deliver high-quality products, precisely meeting and often exceeding customer expectations for safety, performance, and reliability.




EN series



EN-3M





EN-2M



EN-3M

Foldable Electric Home Care Bed



Trend./Rev. Trend. $\pm 11^\circ$



Chair Position



Wooden Head/Foot Board





3/4 Length Siderail, Rotary Armrest Siderail, Plug-in Siderail

EN-2M

Assembly Electric Home Care Bed

5 Cartons Pack



Home Care Bed



Trend./Rev. Trend. $\pm 11^\circ$




Chair Position



Wooden Head/Foot Board



Dust-Proof Lockable Castor



Mattress Stopper


3 Siderails



Rotary Armrest Siderail
 Siderail installs on bedside, caregivers could adjust angles upon patient's need.



Plug-in Siderail
 Siderail installs on bedside or foot-end.



3/4 Length Siderail

Siderails for below models

ES-02JH | EN-3M | EN-2M | ES-18UDS | ES-17UDS



ES-02JH



EN-3M



EN-2M



ES-18UDS



| Specification | EN-3M | EN-2M | Features | EN-3M | EN-2M |
|------------------------------|-----------------------------|-----------------------------|---|-------------------------------------|-------------------------------------|
| Overall bed length | L : 2140mm | L : 2140mm | Chair Position | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> |
| Overall bed width | W : 970mm | W : 1070mm | Foldable Type | <input checked="" type="checkbox"/> | |
| Overall mattress base length | L : 1935mm | L : 1900mm | Assembly Type | | <input checked="" type="checkbox"/> |
| Overall mattress base width | W : 775mm | W : 910mm | Remote Controller | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> |
| Overall bed height | 4"Lockable Caster 300~700mm | 4"Lockable Caster 300~700mm | Rotary Armrest Siderail | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> |
| Head Elevation | 0~70° | 0~75° | 3/4 Length Siderail | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> |
| Knee Elevation | 0~28° | 0~15° | Plug-in Siderail | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> |
| Trend./Rev. Trend. | 11° | 11° | Dust-Proof Lockable Castor | <input type="checkbox"/> | <input checked="" type="checkbox"/> |
| Work Load | 170kg | 173.5kg | Lockable Castor | <input checked="" type="checkbox"/> | |
| Castor Diameter | 4" | 4" | IV Pole Hole | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> |
| Actuator | 4 | 3 | Mattress Stopper | | <input checked="" type="checkbox"/> |
| Volts & Frequency | AC110 ; 50/60Hz | AC110V ; 50/60Hz | Backup Battery (option) 9V Alkaline Battery x3 | <input type="checkbox"/> | <input type="checkbox"/> |
| Actuator/Control Box | IP54 | IP54 | | | |

Standard Option

| | | |
|-------------------------------|-----------------|-----------------|
| Applicable mattress size (mm) | L1980xW820xH150 | L1980xW890xH100 |
|-------------------------------|-----------------|-----------------|

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree



Joson-Care Enterprise Co., Ltd.

No.11, Aly. 77, Ln. 180, Hushan St., Guishan Dist., Taoyuan City 333024, Taiwan

☎ +886-3-329-0925

✉ overseas@joson-care.com.tw

🌐 www.joson-care.com

Taichung Office

1F., No. 280, Chengdu Rd., Xitun Dist.,
Taichung City 407403, Taiwan

TEL : +886-4-2316-0743

Tainan Office

No. 732, Sec. 1, Zhongzheng Rd., Rende Dist.,
Tainan City 717010, Taiwan

TEL : +886-6-266-2809



Document No.: JCQB-031
Release Date: 2024.06.01
Version: L