



# Joson-Care 2024




Safety Quality Sustainability






Business  
Concept




History




Established in 2000, Joson-Care Co., Ltd. is Taiwan's premier manufacturer of hospital beds and ward equipment. We uphold rigorous quality control standards to deliver exceptional accuracy, durability, and reliability, while maximizing the use of renewable resources.





So far, we have sold our products to more than 40 countries. Such as U.S.A, Japan, France, South Korea, India, Iran, Panama, Thailand, Turkey and U.A.E.

Joson-Care's manufacturing system is meticulously designed to adhere to ISO 13485 and Quality Management System (QMS) standards. This rigorous approach ensures that our processes consistently deliver high-quality products, precisely meeting and often exceeding customer expectations for safety, performance, and reliability.



# 02.16 series




## ES-02JH



## ES-16N





## MS-16NS



# ES-02JH

## Electric Nursing Bed

Standards: IEC 60601-1-2, IEC 60601-2-52



Rotary Armrest Siderail (Option)



Locking mechanism for fixed angles.



Detachable Head/Foot Board



Siderail Holder



Remote Controller



IV Pole Hole




Mattress Stopper

# ES-16N series

## Electric Nursing Bed

Standards: IEC 60601-1-2, IEC 60601-2-52

3/4 Length Siderail



Hospital Bed




3/4 Length Siderail + P-Type Siderail



Attendant Control Panel (Option)



Remote Controller



Detachable Head/Foot Board




Mattress Stopper


# MS-16NS

## Manual Hospital Bed


Standards: IEC 60601-1-2, IEC 60601-2-52



3/4 Length Siderail + P-Type Siderail




Manual Screw drive



IV Pole Hole



Detachable Head/Foot Board



Mattress Stopper





Specification	ES-02JH	ES-16N	MS-16NS
Overall bed length	L : 2130mm	L : 2115mm	L : 2115mm
Overall bed width	W : 1000mm	W : 1000mm	W : 1000mm
Overall mattress base length	L : 1910mm	L : 1910mm	L : 1910mm
Overall mattress base width	W : 845mm	W : 845mm	W : 845mm
Overall bed height	H : 320~630mm	H : 350~660mm	H : 350~660mm
Head Elevation	0~75°	0~85°	0~85°
Knee Elevation	0~30°	0~40°	0~40°
Work Load	170kg	225kg	225kg
Castor Diameter	4"	5"	5"
Actuator	3	3	Manual Crank
Volts	AC 120V	AC 120V	None
Frequency	50/60Hz	50/60Hz	None
Backup Battery	Optional	Optional	None
Actuator/Control Box	IP54	IP54	None
Safety Standard	IEC 60601-2-52	IEC 60601-2-52	IEC 60601-2-52

Applicable mattress size (mm)	L1930xW840xH100	L1930xW840xH100	L1930xW840xH100

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree (ES-02JH ±3°)

Features	ES-02JH	ES-16N	MS-16NS
Remote Controller	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Siderail	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Bumper	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Attendant Control Panel	<input type="checkbox"/>	<input type="checkbox"/>	
Lock-Out Function		<input type="checkbox"/>	
Auto Regression		<input type="checkbox"/>	
Pneumatic Cylinder assisted Hi-Lo Mechanism			<input checked="" type="checkbox"/>
Dust-Proof Central Braking /Steering Module	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Dust-Proof Lockable Castor	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Lockable Castor		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
IV Pole Hole	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Drainage Bag Hanger	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Mattress Stopper	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

Standard     Option



### Joson-Care Enterprise Co., Ltd.

No.11, Aly. 77, Ln. 180, Hushan St., Guishan Dist., Taoyuan City 333024, Taiwan

☎ +886-3-329-0925

✉ overseas@joson-care.com.tw

🌐 www.joson-care.com

#### Taichung Office

1F., No. 280, Chengdu Rd., Xitun Dist.,  
Taichung City 407403, Taiwan

TEL : +886-4-2316-0743

#### Tainan Office

No. 732, Sec. 1, Zhongzheng Rd., Rende Dist.,  
Tainan City 717010, Taiwan

TEL : +886-6-266-2809



Document No.: JCQB-031  
Release Date: 2024.06.01  
Version: L