



Joson-Care 2024




Safety Quality Sustainability






Business
Concept




History




Established in 2000, Joson-Care Co., Ltd. is Taiwan's premier manufacturer of hospital beds and ward equipment. We uphold rigorous quality control standards to deliver exceptional accuracy, durability, and reliability, while maximizing the use of renewable resources.





So far, we have sold our products to more than 40 countries. Such as U.S.A, Japan, France, South Korea, India, Iran, Panama, Thailand, Turkey and U.A.E.


Joson-Care's manufacturing system is meticulously designed to adhere to ISO 13485 and Quality Management System (QMS) standards. This rigorous approach ensures that our processes consistently deliver high-quality products, precisely meeting and often exceeding customer expectations for safety, performance, and reliability.




05.99 series




ES-05HD / ES-05HDS



Head Elevation
0 ~ 85°



Knee Elevation
0 ~ 32°



Bed Elevation



Auto-chair




Safety Standard



225KG

Work Load



3PCS

Actuator




Backup Battery



IP66

Actuator/Control
Box


MS-05HS



Head Elevation
0 ~ 85°



Knee Elevation
0 ~ 32°




Bed Elevation



225KG

Work Load

ES-99HDS



Head Elevation
0 ~ 85°



Knee Elevation
0 ~ 32°




Bed Elevation




225KG

Work Load



3PCS

Actuator



Backup Battery



Auto-chair



Safety Standard





IP66

Actuator/Control
Box

ES-05HD

Electric Hospital Bed

Standards: IEC 60601-1-2, IEC 60601-2-38



Attendant Control Panel (ACP)



Swivel Siderail





Backrest Angle Indicator (Option)

ES-05HDS

Electric Hospital Bed

Standards: IEC 60601-1-2, IEC 60601-2-38



Hospital Bed



Siderails up Bed width 100cm fit to small ward.



Detachable Head/Foot Board



IV Pole Hole





Drainage Bag Hanger



Mattress Stopper

MS-05HS



Manual Hospital Bed




ES-99HDS

Electric Hospital Bed

Standard: IEC 60601-1-2, IEC 60601-2-38



Manual Crank Screw drive



Mattress Stopper



Swivel Siderail



IV Pole Hole



Restraint Belt Hole



Specification	ES-05HD	ES-05HDS	MS-05HS	ES-99HDS
Overall bed length	L : 2170mm	L : 2115mm	L : 2170mm	L : 2115mm
Overall bed width (Siderails down/up)	W : 940mm 1090mm	W : 910mm 1000mm	W : 940mm 1090mm	W : 930mm 1090mm
Overall mattress base length	L : 1940mm	L : 1940mm	L : 1940mm	L : 1940mm
Overall mattress base width	W : 925mm	W : 845mm	W : 925mm	W : 925mm
Overall bed height (Central Braking/Normal Castor)	H : 455~755mm H : 395~695mm	H : 455~755mm H : 395~695mm	H : 455~755mm H : 395~695mm	H : 455~755mm H : 395~695mm
Head Elevation	0~85°	0~85°	0~85°	0~85°
Knee Elevation	0~32°	0~32°	0~32°	0~32°
Trend./Rev. Trend.	None	None	None	None
Work Load	225kg	225kg	225kg	225kg
Caster Diameter	5"	5"	5"	5"
Actuator	3	3	Manual Crank	3
Volts & Frequency	AC 100-240V ; 50/60Hz	AC 100-240V ; 50/60Hz	None	AC 100-240V ; 50/60Hz
Backup Battery	DC 24V 2.0Ah	DC 24V 2.0Ah	None	DC 24V 2.0Ah
Actuator/Control Box	IP66	IP66	None	IP66
Safety Standard	IEC 60601-2-38	IEC 60601-2-38	None	IEC 60601-2-38
Applicable mattress size (mm)	L1900xW890xH100	L1930xW840xH100	L1900xW890xH100	L1900xW890xH100

±20mm tolerance in all dimensions ±30mm in siderail width +5° -0° tolerance in angle degree

Features	ES-05HD	ES-05HDS	MS-05HS	ES-99HDS
Outer/Inner Siderail Control	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Backrest, Trend./Rev. Trend. Angle Indicator	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
CPR Quick Release	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Bumper	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Attendant Control Panel	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Electric CPR	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Remote Controller	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Lock-Out Function	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Auto Regression	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pneumatic Cylinder assisted Hi-Lo Mechanism	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Manual Trend./Rev. Trend.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Dust-Proof Central Braking/Steering Module	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Dust-Proof Lockable Castor	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Lockable Castor	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
IV Pole Hole	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Drainage Bag Hanger	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Mattress Stopper	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

Standard Option



Joson-Care Enterprise Co., Ltd.

No.11, Aly. 77, Ln. 180, Hushan St., Guishan Dist., Taoyuan City 333024, Taiwan

☎ +886-3-329-0925

✉ overseas@joson-care.com.tw

🌐 www.joson-care.com

Taichung Office

1F., No. 280, Chengdu Rd., Xitun Dist.,
Taichung City 407403, Taiwan

TEL : +886-4-2316-0743

Tainan Office

No. 732, Sec. 1, Zhongzheng Rd., Rende Dist.,
Tainan City 717010, Taiwan

TEL : +886-6-266-2809



Document No.: JCQB-031
Release Date: 2024.06.01
Version: L